



# Training course in Bioinformatic tools



Instituto de Microbiología  
Centro de Bioinformática

## **Virtual training**

Module 3: SARS-CoV-2 analysis

May 2023


# COVID-19 pandemics




<sup>1</sup> American Society for Microbiology. (2020). *COVID-19 Resources*. <sup>2</sup> WHO. (2021). *COVID-19 Dashboard*. <sup>3</sup> Ritchie et al. (2021). "Coronavirus Pandemic (COVID-19)" - OurWorldInData.org.

# SARS-CoV-2 virus

Family: *Coronaviridae*  
Genre: *Betacoronavirus*  
Subgenre: *Sarbecovirus*



Source: BioRender.com



Source: Gorbalenya et al. (2020). *Nature Microbiology*.

Figure 1 is a scatter plot showing the relationship between genome size and mutation rate per site per replication for various biological entities. The x-axis represents 'Genome Size' on a logarithmic scale from 2 to 10. The y-axis represents 'Mutation rate (per site/replication)' from -11 to 0. The data points are as follows:

Entity	Genome Size (approx.)	Mutation rate (per site/replication) (approx.)
Viroides	2.8	-2.5
ssRNA, ssDNA	3.8	-4.5
dsDNA	4.8	-7.5
Mouse	8.5	-10.5



Source: Duffy *et al.* (2018). *PLOS Biology*; <sup>2</sup> Borges *et al.* (2021).  
*Evolutionary Biology*



# SARS-CoV-2 virus

WHO	PANGO Lineage	RBD	Furin site	Other	Feature
-	B.1	-	-	D614G	↑ transmission
Alpha	B.1.1.7	(E484K), (S494P), N501Y	P681H	D614G, A570D, T716I, S982A, D1118H, (K1991N)	↑ transmission, ↑ viral load, moderate ab evasion
Beta	B.1.351	E484K, K417N, N501Y	-	D614G, A701V	↑ transmission, moderate ab evasion
Gamma	P.1	E484K, K417T, N501Y	-	D614G, H655Y, T1027I	↑ transmission, moderate ab evasion, ↓ vaccine efficiency
Delta	B.1.617.2	(K417N), L452R, T487K	P681R	D614G, D950N	↑ transmission, moderate ab evasion, ↓ vaccine efficiency


Source: Jackson *et al.* (2021). *Nature Reviews*.




Source: Hoffman *et al.* (2020). *Molecular Cell*.



# SARS-CoV-2 lineages and variants



Source: Hodcroft, 2022




10 Apr 2023 - 24 Apr 2023

Variant	Num seq	Freq
23A (Omicron)	466	0.64
22F (Omicron)	170	0.24
22D (Omicron)	56	0.08
23B (Omicron)	17	0.02




# Genomic surveillance

## What is GENOMIC SURVEILLANCE?






# Workflow




# Workflow


## a. RNA extraction



## b. PCR and sequencing library




## c. Sequencing




# RNA to cDNA

---




# Amplicon-based sequencing



# Amplicon-based sequencing

## cDNA amplification (VarSkip Short v2 or NEBNext ARTIC primers)



## Combination of pools




RNA DNA

# SARS-CoV-2 whole-genome primers

New England Biolabs


## UPDATED PRIMERS

### OMICRON VARIANT CLINICAL SAMPLE



# SARS-CoV-2 whole-genome primers


## ARTIC network






# Agarose gel electrophoresis

---




# Workflow


## a. RNA extraction





## b. PCR and sequencing library



## c. Sequencing





# Library preparation - Nanopore




BC01	BC02	BC03	BC04	BC05	BC06	BC07	BC08	BC09	BC10	BC11	BC12	A
BC13	BC14	BC15	BC16	BC17	BC18	BC19	BC20	BC21	BC22	BC23	BC24	B
BC25	BC26	BC27	BC28	BC29	BC30	BC31	BC32	BC33	BC34	BC35	BC36	C
BC37	BC38	BC39	BC40	BC41	BC42	BC43	BC44	BC45	BC46	BC47	BC48	D
BC49	BC50	BC51	BC52	BC53	BC54	BC55	BC56	BC57	BC58	BC59	BC60	E
BC61	BC62	BC63	BC64	BC65	BC66	BC67	BC68	BC69	BC70	BC71	BC72	F
BC73	BC74	BC75	BC76	BC77	BC78	BC79	BC80	BC81	BC82	BC83	BC84	G
BC85	BC86	BC87	BC88	BC89	BC90	BC91	BC92	BC93	BC94	BC95	BC96	H
1	2	3	4	5	6	7	8	9	10	11	12	

BCxx: Barcode xx

# Flowcell loading



# MinKNOW software



# EPI2ME software

EPI2ME Labs

EPI2ME Labs



## Install workflows



### epi2me-labs/wf-alignment

This workflow provides an easy way to align Oxford Nanopore reads and gather mapping stats either locally for small amounts of data or at scale in a distributed environment such as a cluster.

Oxford Nanopore Technologies



### epi2me-labs/wf-metagenomics

Identification of the origin of single reads from both both amplicon-targeted and shotgun metagenomics sequencing.

Oxford Nanopore Technologies



### epi2me-labs/wf-artic

This workflow uses implements the ARTIC SARS-CoV-2 methodology on multiplexed MinION, GridION, and PromethION runs.

Oxford Nanopore Technologies




### epi2me-labs/wf-bacterial-genomes

This workflow can be used for a variety of analyses of bacterial genomes including assembly, variant calling, and annotation.

Oxford Nanopore Technologies



# ARTIC Pipeline





# Workflow

## Lineage assignment

Assignment of clades and lineages



Mutations and substitutions



## Phylogeny

Alignment

MAFFT version 7




- Wuhan-Hu-1 reference genome
- Aliview: check alignment

Phylogenetic tree construction



- Maximum Likelihood, GTR, 1000x
- Visualization in iTOL

## Statistics



# Nextclade - Nextstrain

Nextclade<sub>v2.13.0</sub>

Clade assignment, mutation calling, and sequence quality checks



☰ What's new

◀ Back

Done. Total sequences: 64. Succeeded: 64


ID	Sequence name	QC	Clade	Mut.	non-ACGTN	Ns	Gaps	Nucleotide sequence
?	?	?	?	?	?	?	?	?
0	✓ USA/MA-CDCBI-CRSP_4MU27NEAFJHEV65I/2	N M P C	21F (Iota)	31	0	347	33	
1	✓ USA/MA-CDCBI-CRSP_4UWPG2VONGPMHNX	N M P C	21F (Iota)	32	0	0	10	
2	✓ USA/MA-CDCBI-CRSP_4YUUYAIF675H52XC/2	N M P C	20I (Alpha, V1)	35	0	0	19	
3	✓ USA/MA-CDCBI-CRSP_52WHU54AB4LTIWPX/	N M P C	20I (Alpha, V1)	38	0	0	19	
4	✓ USA/MA-CDCBI-CRSP_5EKJSGVTDNKKI7TN/2	N M P C	20I (Alpha, V1)	35	0	480	19	
5	✓ USA/MA-CDCBI-CRSP_5I6Y4GA3ZXMAC2N6/2	N M P C	21A (Delta)	36	0	246	13	
6	✓ USA/MA-CDCBI-CRSP_5MEECONUJW52G5PA	N M P C	20J (Gamma, V3)	38	0	382	9	
7	✓ USA/MA-CDCBI-CRSP_5YGWYJ5EJQ7H6PNU/	N M P C	20I (Alpha, V1)	38	0	233	19	
8	✓ USA/MA-CDCBI-CRSP_62POGJA32XACPJIP/20	N M P C	21C (Epsilon)	29	0	444	0	
9	✓ USA/MA-CDCBI-CRSP_6B7IU3KSFJ3BVLSL/20	N M P C	20I (Alpha, V1)	36	0	0	19	

# Nextclade




Built with [nextstrain/ncov](#). Maintained by the [Nextstrain team](#). Enabled by data from [GISAID](#).

Showing 2764 of 2764 genomes sampled between Dec 2019 and May 2023.




# Nextclade

Geography



# Nextclade

Frequencies (colored by Clade)



- 20I (Alpha, V1) OR B.1.1.7
- 20H (Beta, V2) OR B.1.351
- 20J (Gamma, V3) OR P.1
- 21A (Delta) OR B.1.617.2
- o Omicron: 21K (Omicron) OR BA.1, 21L (Omicron) OR BA.2, 22A (Omicron) OR BA.4, 22B (Omicron) OR BA.5, 22C (Omicron) OR BA.2.12.1, 22D (Omicron) OR BA.2.75, 22E (Omicron) OR BQ.1, 22F (Omicron) OR XBB, 23A (Omicron) OR XBB.1.5, 23B (Omicron) OR XBB.1.16




# Pangolin COVID-19 Lineage assigner

---

pangolin


Phylogenetic Assignment  
of Named Global Outbreak  
Lineages



Pangolin was developed to implement the dynamic nomenclature of SARS-CoV-2 lineages, known as the Pango nomenclature. It allows a user to assign a SARS-CoV-2 genome sequence the most likely lineage (Pango lineage) to SARS-CoV-2 query sequences.

It is available as a command line tool and a web application. The web application was developed by the Centre for Genomic Pathogen Surveillance. The command line tool is open source software available under the GNU General Public License v3.0.

# Pangolin COVID-19 Lineage assigner






























Retry Failed Sequences


Reset entries

Upload another file

Help


File name	Sequence name	Lineage	Assignment probability
— FAILED (Click warning icon for more info) 2 sequences			
	cluster.fasta	EDB003	
	cluster.fasta	EDB004	
— ANALYSED (Click tick icon for more info) 8 sequences 			
	cluster.fasta	EDB001 B.1.1.65  	1.0
	cluster.fasta	EDB002 B.1.1.65  	1.0
	cluster.fasta	EDB005 B  	1.0
	cluster.fasta	EDB006 B.1.1.65  	1.0
	cluster.fasta	EDB007 B.1.1.65  	1.0
	cluster.fasta	EDB008 B.1.1.65  	1.0
	cluster.fasta	EDB009 B.1.1.65  	1.0
	cluster.fasta	EDB010 B  	1.0

Pangolin (version v2.0.7, lineages version 2020-08-29) is built by Áine, JT, Verity, Emily and Andrew. Web Application by

 Centre for Genomic Pathogen Surveillance



# GISAIID database



© 2008 - 2023 | [Terms of Use](#) | [Privacy Notice](#) | [Contact](#)

You are logged in as **Sully Marquez** - [logout](#)

Registered UsersEpiFlu™EpiCoV™EpiPox™My Profile

EpiCoV™ | Search | Downloads | Upload

Search ▼Reset filters

EPI\_ISL ID

Virus name

EPI\_SET ID

Location

South America / Ecuador / ...

Host

Collection

to

Submission

to

Clade

all

Lineage

Variant

AA Substitutions

Nucl Mutations

☐ Complete

☐ High coverage

☐ Low coverage excluded

☐ With patient status

☐ Collection date complete

☐ Under investigation

Text Search

<input type="checkbox"/>	Virus name	Passage de	Accession ID	Collection da	Submission L		Length	Host	Location	Originating
<input type="checkbox"/>	hCoV-19/Ecuador/GENS BIO-INSPI-1815-1	Original	EPI_ISL_17582816	2023-04-17	2023-04-28	⬇	29,789	Human	South America /	Laborat
<input type="checkbox"/>	hCoV-19/Ecuador/GENS BIO-INSPI-1839-2	Original	EPI_ISL_17582815	2023-04-20	2023-04-28	⬇	29,515	Human	South America /	Laborat
<input type="checkbox"/>	hCoV-19/Ecuador/GENS BIO-INSPI-1839-2	Original	EPI_ISL_17582814	2023-04-19	2023-04-28	⬆	29,841	Human	South America /	Laborat

Total: 9,794 viruses

<< < 1 2 3 4 5 > >>

EPI\_SETSelectAnalysisDownload